

## URETER SLITTING SHEAR; Ø 1.7MM; FL 600MM

5 Fr.; upwards cutting; semi-rigid

single action; insert only



**SKU:** 583-490-16A

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscapy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Uretero-Renoscapy</a>
<b>TYPE</b>	<a href="#">Semi-rigid Instruments</a>