

URETER SLITTING SHEAR; Ø 1.7MM; FL 600MM

5 Fr.; upwards cutting; semi-rigid

single action; insert only



SKU: 583-490-16A

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments