

URETER SLITTING SHEAR; Ø 1.7MM; FL 600MM

5 Fr.; upwards cutting; semi-rigid

single action; with LuerLock;insert only



with Luer Lock

SKU: 583-490-16LLA

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Uretero-Renoscapy](#)

TYPE [Semi-rigid Instruments](#)