

URETERO- RENOSCOPE

6.0/7.5Fr.; working channel 4Fr.

430mm

SKU: 601-110-20

Categories: [Endoscopy](#), [Equipment](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)