

VAPORIZATION ELECTRODE

Ø 3mm; toothed; 24 / 26Fr.

yellow/blue; bipolar system



SKU: 630-240-45B

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|----------------------------------|
| UNIT | pcs. |
| RANGE | Resectoscopy 4.0 |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|----------------------------------|
| UNIT | pcs. |
| RANGE | Resectoscopy 4.0 |
| TYPE | Electrodes |