

## VAPORIZATION ELECTRODE

Ø 5mm; smooth; 24 / 26Fr.

yellow; monopolar



**SKU:** 630-240-60

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Electrodes</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Electrodes</a>