

VAPORIZATION ELECTRODE

Ø 5mm; smooth; 24 / 26Fr.

yellow/blue; bipolar system



SKU: 630-240-60B

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Electrodes