

## VERES INSUFFLATION CANNULA; HIGH FLOW

Ø 2.7mm; FL 150mm

stopcock; Luer-Lock

**SKU:** 618-515-15

**Categories:** [Endoscopy](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [VERRES Cannulas](#)