

VERES INSUFFLATION CANNULA; HIGH FLOW

Ø 2.7mm; FL 120mm

stopcock; Luer-Lock; stainless steel



SKU: 618-555-12

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [VERRES Cannulas](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [VERRES Cannulas](#)