

WAGENER EAR HOOK

3.5mm

14.0cm (5.6")

SKU: 111-588-02

Categories: [Open Surgery](#), [Otology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
Aesculap	OF283
KLS Martin	36-585-02
AUTHORS	WAGENER
Medicon	62.16.04



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
Aesculap	OF283
KLS Martin	36-585-02
AUTHORS	WAGENER

Medicon

62.16.04