

WAGENER EAR HOOK

1.5mm

14.0cm (5.6")

SKU: 111-588-06

Categories: [Open Surgery](#), [Otology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
KLS Martin	36-585-07
AUTHORS	WAGENER



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Loops, Hooklets
KLS Martin	36-585-07
AUTHORS	WAGENER