

## WAGENER EAR HOOK

4.0mm

14.0cm (5.6")

**SKU:** 111-588-01

**Categories:** [Open Surgery](#), [Otology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Loops, Hooklets</a>
<b>Aesculap</b>	OF282
<b>KLS Martin</b>	36-585-01
<b>AUTHORS</b>	<a href="#">WAGENER</a>
<b>Medicon</b>	62.16.05



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Loops, Hooklets</a>
<b>Aesculap</b>	OF282
<b>KLS Martin</b>	36-585-01
<b>AUTHORS</b>	<a href="#">WAGENER</a>

**Medicon**

62.16.05