

WATER JET PISTOL

complete with 9 tops 250-501-01 - 09

1m pressure tube; strap and suspen. hood



SKU: 250-501-00

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Water jet Pistol](#)

KLS Martin 73-101-00

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Water jet Pistol](#)

KLS Martin 73-101-00