

WAUGH DRESSING FORCEPS

straight; serrated

20.0cm (8"); width of jaw 2.0 +/-0.2

SKU: 020-130-20

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-171-20
AUTHORS	RAMSEY
Medicon	06.42.20
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
KLS Martin	12-171-20
AUTHORS	RAMSEY
Medicon	06.42.20

SPECIFICATIONS

[serrated](#), [straight](#)