

## WEBSTER NEEDLE HOLDER

straight; smooth

13.0cm (5.2")



**SKU:** 050-028-13

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM202
<b>KLS Martin</b>	20-334-13
<b>AUTHORS</b>	<a href="#">WEBSTER</a>
<b>Medicon</b>	10.10.12
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM202
<b>KLS Martin</b>	20-334-13

**AUTHORS**

[WEBSTER](#)

**Medicon**

10.10.12

**SPECIFICATIONS**

[straight](#)