

WEITLANER SELFRETAINING RETRACTOR

blunt; 2X3 teeth

10.5cm (4.2")




SKU: 070-671-10

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV073
KLS Martin	15-749-10
AUTHORS	WEITLANER
Medicon	24.22.11
SPECIFICATIONS	blunt , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV073
KLS Martin	15-749-10

AUTHORS

[WEITLANER](#)

Medicon

24.22.11

SPECIFICATIONS

[blunt, toothed](#)