

WEITLANER SELFRETAINING RETRACTOR

sharp; 2X3 teeth

10.5cm (4.2")




SKU: 070-670-10

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV070
KLS Martin	15-748-10
AUTHORS	WEITLANER
Medicon	24.22.10
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV070
KLS Martin	15-748-10

AUTHORS

[WEITLANER](#)

Medicon

24.22.10

SPECIFICATIONS

[sharp, toothed](#)