

WILDCAT GUIDE WIRE

Platin tip J-form , teflon coated

Ø 0.035", length 650cm single use



SKU: PK-3LQPK0850650XJ-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopia](#)