

## WILDCAT GUIDE WIRE

Platin tip J-form , teflon coated

Ø 0.035", length 650cm single use



**SKU:** PK-3LQPK0850650XJ-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopia](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopia](#)