

WILDCAT GUIDE WIRE

Platin tip, teflon coated single use

Ø 0.035", length 650cm



SKU: PK-3LQPK0850650X-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)