

## WIRE CUTTER FRONT ACTION

soft wire up to Ø3.2mm; hard to Ø2.5mm

26.0cm (10.4")

**SKU:** 190-776-26

**Categories:** [Open Surgery](#), [Orthopaedics](#)



## ZUSÄTZLICHE INFORMATIONEN

|                   |                              |
|-------------------|------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>         |
| <b>TYPE</b>       | <a href="#">Wire Forceps</a> |
| <b>KLS Martin</b> | 22-528-03                    |

## ZUSÄTZLICHE INFORMATIONEN

|                   |                              |
|-------------------|------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>         |
| <b>TYPE</b>       | <a href="#">Wire Forceps</a> |
| <b>KLS Martin</b> | 22-528-03                    |