

## WIRE LOOP FOR TONSIL SNARE

Ø 40 x 0.4mm; package of 100 pieces



**SKU:** 132-520-40

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)


### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Tonsil Snares](#)

**Aesculap** OM757

**KLS Martin** 39-511-40



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Tonsil Snares](#)

**Aesculap** OM757

**KLS Martin** 39-511-40