

WIRE LOOP FOR TONSIL SNARE

Ø 40 x 0.4mm; package of 100 pieces



SKU: 132-520-40

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)


ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tonsil Snares](#)

Aesculap OM757

KLS Martin 39-511-40



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Tonsil Snares](#)

Aesculap OM757

KLS Martin 39-511-40