

WORKING ELEMENT

passive cutting; closed handle

for Ø 2.9mm telescopes

Only for monopolar use



SKU: 630-335-60

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 2.9
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 2.9
TYPE	Working Elements