

WORKING ELEMENT

active cutting; closed handle

for Ø 4mm telescopes;HF-connection above

Only for monopolar use



SKU: 630-235-11

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements