

## WORKING ELEMENT

passive cutting; closed handle

for Ø 4mm telescopes;HF-connection above

Only for monopolar use



**SKU:** 630-235-61

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>