

WORKING ELEMENT

passive cutting; closed handle

for Ø 4mm telescopes;HF-connection above

Only for monopolar use



SKU: 630-235-61

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Resectoscopy 4.0](#)
TYPE [Working Elements](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Resectoscopy 4.0](#)
TYPE [Working Elements](#)