

WORKING ELEMENT; BIPOLAR; TITANIUM BLUE

KS compatible incl. 2 ea cone sleeves

for 4mm/302mm telescopes



SKU: 630-350-05

Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Working Elements