


WORKING ELEMENT; BIPOLAR

passive cutting, for 15.5CH Urethrorome

for Ø2.9mm telescope, HF-connection up



SKU: 630-405-60BU

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Urethrotomy
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Urethrotomy
TYPE	Working Elements