

## WORKING ELEMENT; BIPOLAR

passive cutting, for 15.5CH Urethrorome

for Ø2.9mm telescope, HF-connection up



**SKU:** 630-405-60BU

**Categories:** [Endoscopy](#), [Urology](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Urethrotomy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Urethrotomy</a>
<b>TYPE</b>	<a href="#">Working Elements</a>