

WORKING ELEMENT; PASSIVE; BIPOLAR

18.5 Fr. for Ø 2.9mm telescope; 30°

QL connection; colour code green

for bipolar use



SKU: 582-180-00B

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Ellipse , Mini Resectoscopy (Midi) , Original
TYPE	Working Elements

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Ellipse , Mini Resectoscopy (Midi) , Original
TYPE	Working Elements