


## WORKING ELEMENT PASSIVE TITANIUM BLUE

KS compatible incl. 2 Cone sleeves

for 4mm/302mm Telescope

GYRUS PK



**SKU:** 630-350-55

**Categories:** [Endoscopy](#), [Gynecology](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Resectoscopy 4.0</a>
<b>TYPE</b>	<a href="#">Working Elements</a>