

## WORKING SHEATH WITH 1 INSTRUMENT CHANNEL

5Fr.; for Ø 2.9mm telescopes/0°/300mm

for flexible instruments

for visual urethrotome 15.5Fr. (630-401-15) white



**SKU:** 630-401-30

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Urethrotomy](#)

**TYPE** [Sheaths](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Urethrotomy](#)

**TYPE** [Sheaths](#)