

## WORKING SHEATH WITH

1 working channel 5 Fr.

Quick Lock System for endoscope 2.9mm

0° / 30°

**SKU:** 582-150-39GU

**Categories:** [Endoscopy](#), [GUBBINI/Gyn](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Original</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Original</a>
<b>TYPE</b>	<a href="#">Sheaths</a>