

YANKAUER SUCTION TUBE

removable suction tip; Luer connector

inner Ø 2.0mm; outer Ø 6.0mm; 29.0cm

SKU: 085-000-27

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF944
KLS Martin	18-543-27
AUTHORS	YASARGIL
Medicon	09.85.01

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF944
KLS Martin	18-543-27
AUTHORS	YASARGIL

Medicon

09.85.01