

## YASARGIL APPLYING FORCEPS

straight; for Mini vessel clips

18.0cm (7.2")

**SKU:** 244-610-18

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE548
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.96
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE548
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.96

## SPECIFICATIONS

[jaw rigid](#)