

## YASARGIL APPLYING FORCEPS

bayonet; for Mini vessel clips

21.0cm (8.4")

**SKU:** 244-611-21

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE572
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.98
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE572
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.98

## SPECIFICATIONS

[jaw rigid](#)