

## YASARGIL APPLYING FORCEPS

bayonet; for Mini vessel clips

23.0cm (9.2")

**SKU:** 244-611-23

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE574
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.99
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>Aesculap</b>	FE574
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.56.99

## SPECIFICATIONS

[jaw rigid](#)