

YASARGIL APPLYING FORCEPS; JAW ROTATEABLE

bayonet; for Mini vessel clips

21.0cm (8 1/4")

SKU: 244-613-21

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable