

## YASARGIL APPLYING FORCEPS; JAW ROTATEABLE

bayonet; for Mini vessel clips

23.0cm (9")

**SKU:** 244-613-23

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">jaw rotatable</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">jaw rotatable</a>