

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 9.0cm (3.6")rotatable

for Mini vessel clips



SKU: 244-680-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.92
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.92
SPECIFICATIONS	jaw rigid , shaft rotatable