

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 9.0cm (3.6")rotatable

for Mini vessel clips; TITAN



SKU: 244-680-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | jaw rigid , shaft rotatable |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | jaw rigid , shaft rotatable |