

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 11.0cm (4.4")rotatable

for Mini vessel clips; TITAN



SKU: 244-680-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable