

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 13.0cm (5.2")rotatable

for Mini vessel clips



SKU: 244-680-13

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| AUTHORS | YASARGIL |
| Medicon | 58.37.94 |
| SPECIFICATIONS | jaw rigid , shaft rotatable |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| AUTHORS | YASARGIL |
| Medicon | 58.37.94 |
| SPECIFICATIONS | jaw rigid , shaft rotatable |