

## YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 13.0cm (5.2")rotatable

for Mini vessel clips; TITAN



**SKU:** 244-680-13T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a> , <a href="#">shaft rotatable</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a> , <a href="#">shaft rotatable</a>