

## YASARGIL APPLYING FORCEPS; RING HANDLE

jaw and shaft 7.0cm (2.8")rotatable

for Mini vessel clips



**SKU:** 244-682-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Applying Forceps for Aneurysm Clips</a>             |
| <b>AUTHORS</b>        | <a href="#">YASARGIL</a>  |
| <b>Medicon</b>        | 58.34.91  |
| <b>SPECIFICATIONS</b> | <a href="#">jaw rotatable</a> , <a href="#">shaft rotatable</a> |

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Applying Forceps for Aneurysm Clips</a>             |
| <b>AUTHORS</b>        | <a href="#">YASARGIL</a>  |
| <b>Medicon</b>        | 58.34.91  |
| <b>SPECIFICATIONS</b> | <a href="#">jaw rotatable</a> , <a href="#">shaft rotatable</a> |