

## YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 11.0cm (4.4")rotatable

for Standard vessel clips



**SKU:** 244-640-11

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.32.93
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a> , <a href="#">shaft rotatable</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Applying Forceps for Aneurysm Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.32.93
<b>SPECIFICATIONS</b>	<a href="#">jaw rigid</a> , <a href="#">shaft rotatable</a>