

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 7.0cm (2.8")rotatable

for Mini vessel clips

SKU: 244-680-07

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.91
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
AUTHORS	YASARGIL
Medicon	58.37.91
SPECIFICATIONS	jaw rigid , shaft rotatable