

YASARGIL APPLYING FORCEPS; ROUND HANDLE

rotatable jaw; shaft 11.0cm (4.4")

for Standard vessel clips; TITAN

SKU: 244-637-11T

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | jaw rotatable |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Applying Forceps for Aneurysm Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| SPECIFICATIONS | jaw rotatable |