


YASARGIL MINI VESSEL CLIP

permanent; 7mm; opening 4.4mm

110 gms; curved



SKU: 244-467-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

Aesculap FE722

AUTHORS [YASARGIL](#)

Medicon 58.57.07

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE722

AUTHORS [YASARGIL](#)

Medicon 58.57.07

SPECIFICATIONS

[curved](#), [permanent](#)